

Overcoming technical difficulties of *in-vitro* assay-development.




Fiona Leslie, Chloe Whiting, Megan Webster




Newcells Biotech Limited, The Biosphere, Draymans Way, Newcastle Helix, Newcastle upon Tyne. UK NE4 5BX

Modelling Fibroblast-to-Myofibroblast Transition (FMT)

- TGF- β 1 mediated activation of fibroblasts, their transition to α -SMA expressing myofibroblasts, and associated increase in the expression and deposition of extracellular matrix proteins represents one of the underlying pathologic mechanisms of fibrosis.
- Modelling this process *in-vitro* allows the testing of potential anti-fibrotic therapeutics for their ability to reduce or reverse FMT.
- The development and validation of robust assay systems requires significant time and effort. Our team worked to overcome multiple technical challenges to optimise our high-throughput, high-content imaging FMT-assay. Here we highlight some key difficulties and resolved outcomes.



Optimization of Plate Coatings



To ensure consistent adherence of cells to culture wells throughout the FMT-assay process, including mechanical washing, we optimized cell seeding densities and culture media to control for cell proliferation (data not shown). After observing significant cell detachment due to mechanical washing, we compared plate coatings to promote fibroblast attachment and reduce assay variation. Further work was performed to optimized immunocytochemical detection of both α SMA and collagen I.



Optimization of Plate Coatings

Collagen I	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	
A	4633	360102	376803	357530	297886	317401	325921	359880	305262	364636	310480	379652	174563	330468	229237	384371	281691	337497	316636	364475	359855	330675	367648		
B	16	272484	34410	283351	284245	351071	351950	420918	489174	635092	767552	865093	3444043	1130072	1433263	1051188	1185452	1249012	1143987	253380	330708	353473	338874		
C	8385	520321	42303	495519	379763	466734	466664	466664	466664	466664	466664	466664	466664	466664	466664	466664	466664	466664	466664	466664	466664	466664	466664		
D	0	469454	459567	436511	245674	381946	338427	462147	370836	751738	784975	1026314	1098517	1318065	1232523	1083192	1151910	1192048	1212304	276215	408110	338488	357207		
E	21	853607	707854	941663	282607	366240	321529	382983	363348	507867	551417	1089911	1472189	1639538	1159491	1119638	1082603	1109849	1092655	263415	543768	65845	563684		
F	0	1108002	62440	843896	258357	275207	735287	735287	735287	735287	735287	735287	735287	735287	735287	735287	735287	735287	735287	735287	735287	735287	735287		
G	0	1226594	82386	1010101	343001	321771	344687	331509	350282	515995	559004	965888	1277134	1055511	1064798	1057162	1098613	1024389	1065289	294162	789513	86712	957399		
H	0	29842	1588992	62895	914375	264243	379832	359033	358451	397561	692818	791872	636833	856870	915768	950782	893949	950396	1071572	1039137	232072	496246	31149	839994	
I	191409	1179077	116336	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706		
J	222188	1260300	103950	1200143	378853	365778	350147	414549	387146	770361	693754	804532	797821	1012703	971311	1059127	951964	883631	964787	228044	1036274	102333	907661		
K	177907	1246709	106392	1089907	308139	321324	343088	334730	303012	577052	519399	891798	1123668	1270865	1319397	1111476	1119724	997700	1018140	283169	1140993	105873	1141343		
L	0	1181769	93897	1094163	264428	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706		
M	0	1181769	93897	1094163	264428	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706		
N	0	1309429	33364	852190	299820	240846	259964	244450	341527	534334	526545	783785	1317273	1328636	1187583	1235584	971261	1025072	852505	288026	988804	84886	1011153		
O	0	1078883	63375	1006435	328528	231208	257474	223995	315550	330499	334706	269599	258442	346454	319082	365100	217365	254486	274515	220332	1173602	101301	1233786		

α SMA	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A	2064	11253	911	7165	20371	20269	13838	9878	18311	21862	13685	22088	44679	25252	57044	14652	5208	12211	9976	13297	17425	12263	15028	
B	1153	8782	11368	25716	19882	14046	13937	9193	13355	43132	56915	133032	59892	491431	152340	710302	562485	578133	622502	8055	14911	19145	8080	
C	271	8690	30392	8940	18978	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	
D	15575	33597	14002	7506	52396	5566	13207	19974	15283	34728	24258	98170	80445	348444	202411	444765	664016	591296	603184	10600	11210	13141	23614	
E	22135	58189	6288	21153	18003	11643	19512	19229	18170	24377	18514	84562	55018	371784	237546	598096	597534	335448	604103	11119	17495	13128	13711	
F	16601	35543	7674	24186	32518	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	
G	8293	50744	10727	47477	37414	9550	20506	7364	24023	24737	29233	59133	70919	308889	403383	454565	547393	416167	584612	12812	77873	1025	44996	
H	2852	60169	5292	63950	28155	16612	13916	14434	27482	41561	47228	172091	127149	407085	506347	584320	656136	620251	620065	31938	66584	2461	68952	
I	980	358397	46441	304646	24288	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	
J	1903	98819	40704	493849	15536	21734	15578	15348	18554	57481	30882	142724	146274	396527	348814	665861	581931	528135	500484	17307	41343	13764	191593	
K	3760	116605	80940	577787	17315	5316	20667	7258	19919	21917	22086	126884	82987	479135	239554	676135	627486	567126	600388	18768	499991	10707	168564	
L	678	265238	39116	424374	26768	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	1169204	2561706	
M	5293	109821	30943	389674	23842	22749	16510	19254	21042	34662	24160	79104	75000	270841	185517	620401	574862	629898	479372	4557	509399	16104	232898	
N	721	422566	634580	381880	10567	26299	14628	18025	21342	5640	20804	9140	15370	6194	26091	4036	7182	4556	7194	9896	482858	26754	369576	

To ensure consistent data, independent of well position on the plate, α SMA and collagen I were quantified following the stimulation of primary fibroblasts with increasing TGF- β 1 concentration across both the horizontal and vertical axis (arrows). The effects of two different stimulation methods were also tested (dashes). Applied conditional formatting shows increased expression of α SMA and collagen I, ranging green through red respectively.

Optimization of Culture Media



Ensuring a consistent, substantial, assay window is vital to establishing a reproducible *in-vitro* assay for compound screening. To promote the activation of primary fibroblasts following stimulation with TGF- β 1 we worked to optimize assay culture media. As shown above, media 3 promoted the expression of both α SMA and collagen1 at the gene and protein level. Together with our refined segmental analysis masks we can accurately detect changes in the expression and deposition of α SMA and collagen I in response to treatment with potential therapeutic compounds.

Newcells Biotech's FMT-assay

- Newcells' FMT-assay utilises high-content imaging to determine the effects of potential therapeutics on fibroblast activation and collagen expression and deposition in high-throughput format.

Assay Features:

- 3 validated primary HLF donors
- High-throughput 384-well format
- 6-point dose response for up to 6 TAs per plate
- 6 technical replicates per condition
- Plus validated experimental assay controls

Available Assay Readouts:

- ✓ Cell Number
- ✓ Collagen I
- ✓ α SMA

Further details available online or contact us at enquiries@newcellsbiotech.co.uk



Also available: Small Airway Epithelial Cell (SAEC) model.

